SENSOR RESERVOIR® / SENSOR VORKAMMER®

Ⓛ Patientenhandbuch | ⓊⓀⓁ USA Patient Manual

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THE COMPANY

Christoph Miethke GmbH & Co. KG is a company based in Berlin-Brandenburg region that develops, manufactures and markets innovative neurosurgical implants for the treatment of hydrocephalus. In the course of our work we established successful partnerships with numerous hospitals worldwide. The purpose of this booklet is to provide you and your family with some understanding of the treatment of hydrocephalus.

The successful treatment of this condition has only been possible since the 1950s. In a dramatic race against time to save the life of his son, Casey, who suffered from hydrocephalus, a technician named John D. Holter developed, in only a few weeks, a novel silicone valve. Despite the fact that, since its first implantation in March 1956, this valve has proven to be clinically effective and a giant step in the treatment of this condition, there are many patients today who experience considerable problems with hydrocephalus valve systems.

Christoph Miethke utilised the knowledge gained in 50 years of valve treatment and developed a new generation of highprecision valves made of the metal titanium. For the first time, there are valve systems available that consistently take into account the physical conditions of brain fluid drainage and can thus maintain a physiological brain pressure, independent of the body position of the patient.

Fig. 1: Anatomic sketch of the cranium
(see inner cover page)

- 1) skull
- 2) brain
- 3) cerebrospinal fluid
- 4) lateral ventricle
- 5) third ventricle
- 6) fourth ventricle

BASIC ANATOMY

The human brain (fig. 1) is surrounded by a special substance known as cerebrospinal fluid (CSF). Cerebrospinal fluid is produced in several chambers, so-called ventricles, that are found within the brain. The channels, by which the ventricles are interconnected, constitute a complex drainage system. The fluid in the brain circulates through these ventricles and eventually flows into the venous blood. The function of this fluid is to protect the brain from mechanical damage. The CSF also regulates the internal brain pressure (intracranial pressure, ICP), keeps the brain tissue moist and transports the products of metabolism.

CLINICAL PICTURE OF THE CONDITION

In healthy people, a balance exists between the production and resorption of cerebrospinal fluid. In infants, approx. 100 ml of this fluid is produced every day; in small children, the daily production is approx. 250 ml, in grown-ups approx. 500 ml. If the amount of fluid produced exceeds the amount resorbed, the ventricles expand, leading to the condition known as hydrocephalus (fig. 2). The term hydrocephalus refers to the continuous increase of the volume of "water" (hydro) in the "head" (cephalus). This condition is often observed at birth (congenital hydrocephalus), but it can also develop later in life, e.g., as the result of inflammation, hemorrhage or severe head injury, or after brain surgery. Such cases are referred to as acquired hydrocephalus.


Fig. 2: Ventricle size
a) normal, b) hydrocephalus

A further distinction is made between obstructive hydrocephalus and communicating hydrocephalus. In obstructive hydrocephalus, the links between the ventricles of the brain are interrupted so that the ventricles cannot "communicate" with each other. Cases in which the ventricles are interlinked through open channels, but resorption of cerebrospinal fluid is impaired, are diagnosed as communicating hydrocephalus.

CLINICAL SYMPTOMS OF THE CONDITION

In infants, the cranial bones have not yet grown together solidly. The increasing volume of cerebrospinal fluid causes the head to increase in circumference while, at the same time, brain tissue disintegrates. From the age of about 2, the hardened skull prevents any growth of the head's circumference. In this case, the increase in fluid volume leads to a massive pressure increase, resulting in the expansion of the brain ventricles and the com-

pression of the brain itself. The consequence for infants and adults can be irreversible brain damage. Symptoms (depending on the severity of the disorder) include nausea, headache, vomiting, impaired coordination, drowsiness and, in the end, unconsciousness.

METHODS OF TREATMENT

An alternative in the treatment of hydrocephalus is the implantation of a shunt system which guarantees the drainage of CSF from the ventricles to the peritoneum. Such a shunt system consists of several catheters and a hydrocephalus valve.

Additionally, a *SENSOR RESERVOIR* and a *SENSOR PRECHAMBER* can be integrated, thereby fulfilling the function of a conventional reservoir:

1. The reservoir allows the measurement of intraventricular pressure by puncture, the injection of medicine and the control of the valve's functionality. The robust titanium housing of the integrated measuring cell prevents a perforation of the reservoir.

2. The *SENSOR RESERVOIR* and the *SENSOR PRECHAMBER* (fig.1.) offer a non-invasive possibility to carry out a functionality control of the Shunt(-system). Therefore it might be of particular interest to understand the pressure setting inside the reservoir and the change in the pressure settings. Due to the integrated measuring cell inside the reservoir, the *SENSOR RESERVOIR* and the *SENSOR PRECHAMBER* enable the measurement of these pressures, to record these results and to consequently draw conclusions from this data, which helps to access the functional efficiency of the shunt system.

The *SENSOR RESERVOIR* and the *SENSOR PRECHAMBER*, as well as the catheters, consist of the non-magnetic materials PEEK, titanium and silicone.

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt.

MRI examinations can produce artifacts.

Should the patient be carrying a Sensor Reservoir a MRT-investigation is possible without damage to the device and patient. The *SENSOR RESERVOIR*, the *SENSOR PRECHAMBER* and all Miethke Shunt systems are designed in such a way that they are able to withstand magnetic energies up to 3 Tesla.


Fig.3: Functionality of

- a) the *SENSOR RESERVOIR* and
b) the *SENSOR PRECHAMBER* with integrated measuring cell

OPERATIONAL FUNCTIONALITY OF THE *SENSOR RESERVOIRS/ SENSOR PRECHAMBER*

The data of the measuring cell is scanned and then displayed by using the Reader Unit of the *SENSOR RESERVOIR/SENSOR PRECHAMBER*. The SD-card automatically stores the measurement data and a later data analysis is then made possible. The functionality of the shunt can be determined by the relative curve of the CSF pressure inside the shunt. Therefore occlusions of the shunt can be detected in a non-invasive way and even the

mechanical functional loss of the shunt can be detected, localized and evaluated.

Fig.4: Pressure measurement with the Reader Unit of the *SENSOR RESERVOIR/ SENSOR PRECHAMBER*

THERAPY COMPLICATIONS

In certain instances, complications can result from the treatment of hydrocephalus with a shunt system. As with any surgical procedure, an infection can develop. Unfortunately, problems can also arise that are either directly or indirectly linked with the implanted valve system. More information can be found in each relevant IFU.

Depending on the *SENSOR RESERVOIR* and the *SENSOR PRECHAMBER* note that frequent pumping can lead to excessive drainage and therefore to unphysiological pressure values. The patient should be informed about the risk. If the patient suffers from an elevated body temperature, an impairment of the Reader Unit is possible. (see IFU Reader Unit for *SENSOR RESERVOIR*). While scanning, a rise of temperature inside the *SENSOR RESERVOIR* or the *SENSOR PRECHAMBER* could occur. In this case the measurement will be stopped at 39°C by a temperature sensor.

Hydrocephalus patients who experience headache, dizziness, unnatural gait or similar symptoms should consult a physician without delay.

WARNING NOTES

The use of energy-emitting devices such as defibrillators and HF devices (high-frequency equipment/ devices) can lead to failure of the measuring cell!

The measuring cell could be damaged in the event that an electrical current from an external source is directed through the body,

Radiotherapy and radionuclide examinations of the patient could lead to a failure of the measuring cell.

Diagnostic ultrasound examinations are safe. Ultrasound therapy could lead to an unintended concentration of the ultrasonic field.

AFTER THE OPERATION

The everyday activities of patients with *SENSOR RESERVOIR* and the *SENSOR PRECHAMBER* are not restricted. However, patients should consult their attending physician before major physical exertion (e. g. hard physical work, strenuous sports). Apart from that, we recommend medical check-ups at regular intervals. If the patient undergoes activities with high pressure and impact load (f.ex: diving, boxing, soccer, etc.), a *SENSOR RESERVOIR* or a *SENSOR PRECHAMBER* should not be used in their shunts.

The patient should avoid knocks or pressure.

PATIENT ID AND SD-CARD

Each *SENSOR RESERVOIR* or the *SENSOR PRECHAMBER* is supplied with an individual patient ID. The ID is filled out by the attending physician and contains important information for the follow-up examinations.

In addition to the patient pass, the patient receives a SD-card. The SD-card contains all individual data concerning the implant. In case the SD-card is lost it can be reordered by specifying the manufacturer's ASIC-ID or serial number of the SD-card.

A BRIEF PATIENT GLOSSARY

Anatomy

A guide to the structure of body components

Arachnoid

Connective tissue in the brain that lies between the dura mater and the pia mater

Catheter

Tube

Cerebrospinal fluid

Watery spinal fluid in the brain

Communicating vessels

Vessels that are connected by a channel

Computed tomography (CT)

Imaging technique whereby „slices“ of the body are recorded with an X-ray scanner

Drainage

Drainage of accumulated fluid

Dura mater

The hardest component of the meninx

Fluid component

Cerebrospinal fluid component

Fontanel

Connective tissue opening in a young infant's skull that later ossifies

Implant

Substance that is placed in the human body to replace a particular function for a limited period of time or for the rest of the patient's life

Lumbar puncture

Puncture of the spinal channel at the lower spine

Lumbo-peritoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain, by way of the region of the lumbar vertebrae in the abdominal cavity

Meninges

Membrane found in the brain and spine

Meningitis

Inflammation of the meninx

Minimally invasive

Minimally infiltrating

Overdrainage

Undesirable outward flow of cerebrospinal fluid

Peritoneum

Membrane that covers the pelvic and abdominal cavities

Pia mater

Component of the soft meninx containing blood vessels

Piaarachnoid

Soft component of the meninx that is divided into the arachnoidea and pia mater

Puncture

Insertion of a hollow needle or a trocar into a vessel for the purpose of removing fluid

Resorption

Suctioning or removal of material through skin, mucosa, or tissue

Shunt

A passage between two channels – here a catheter drainage system with an integrated valve

Spinal column

Element of the central nervous system located within the vertebral channel

Subcutaneous pressure

Pressure beneath the skin

Subdural hematoma

An accumulation of blood between the brain and cranium

Ventricle of the brain

Intracranial space containing cerebrospinal fluid

Ventricular peritoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain directly into the abdominal cavity (abdominal catheter)

FOLLOW-UP EXAMINATIONS

A follow-up examination must be carried out in all cases.

Date	Treatment

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany